

SPX-29

SCALERS, PERIODONTOLOGY

**Efficiency. Precision**  
Care – gentle removal,  
superior control.

*Surgiplex. The Difference You Hold*



[www.surgiplexsrl.com](http://www.surgiplexsrl.com)

ABOUT US

# The Bridge Between Manufacturing Excellence and Global Supply

**S**urgiplex SRL is a Romania-based export and distribution company specializing in the global supply of high-quality surgical instruments. Operating from Bucharest, Surgiplex acts as the commercial and logistics hub of the Surgilabs Group, a partnership designed to unify manufacturing expertise with European-level quality assurance and distribution.

In collaboration with Surgilabs (Pakistan) – our dedicated manufacturing partner – we ensure a seamless connection between production precision and international healthcare delivery. Surgilabs is responsible for crafting each surgical instrument to global standards, while Surgiplex manages export, compliance, and customer relations, ensuring every client receives products that reflect reliability, safety, and performance.

**OUR COLLABORATION  
TRANSFORMS PRECISION  
INTO RELIABILITY – EV-  
ERY INSTRUMENT WE EX-  
PORT CARRIES A PROMISE  
OF PERFORMANCE.**

Together, Surgiplex and Surgilabs form a vertically integrated network – one that transforms local craftsmanship into a global benchmark of quality and care.



To be a trusted global partner in delivering innovative, precision-crafted surgical instruments that improve healthcare outcomes worldwide.

To connect healthcare providers with superior-quality surgical products through ethical practices, seamless logistics, and exceptional customer care, powered by our dedicated production partner, Surgilabs.

**- INTEGRITY**

Transparent and ethical in every process.

**- EXCELLENCE**

Instruments that exceed international standards.

**- INNOVATION**

Modern designs and flexible solutions.

**- ACCESSIBILITY**

Affordable instruments for healthcare worldwide.

**- PARTNERSHIP**

Building long-term networks with distributors and providers.

20+

Active distributors and partners worldwide.

5-7%

Annual revenue investment on R&D

Precision. Partnership. Progress.

PRODUCT RANGE

# From Design to Delivery — A Complete Range of Surgical Solutions

At Surgiplex SRL, our product range reflects a commitment to quality, innovation, and precision across every discipline of surgery.

In partnership with Surgilabs (Pakistan), we manufacture and supply a comprehensive line of surgical, dental, and specialty instruments, engineered for performance, reliability, and ergonomics.

■ **Precision:**

Every edge, hinge, and grip is engineered to meet surgical accuracy.

■ **Performance:**

Designed for longevity, balance, and ease of handling.

■ **Safety:**

Manufactured under strict hygiene and sterilization protocols.

■ **Innovation:**

Continuous design refinement based on global medical feedback.



GENERAL SURGERY INSTRUMENTS  
DENTAL INSTRUMENTS  
ORTHOPEDIC INSTRUMENTS  
GYNECOLOGY & OBSTETRICS  
ENT & OPHTHALMIC INSTRUMENTS  
SINGLE-USE INSTRUMENTS  
CUSTOM & OEM INSTRUMENTS

100+ Total Categories

- **General Surgery Instruments**  
Dissecting & Dressing Forceps  
Surgical Scissors & Shears  
Needle Holders & Clamps  
Retractors (Self-Retaining & Hand-Held)  
Hemostatic & Artery Forceps  
Towel Clips & Sponge Forceps  
Operating & Tissue Forceps  
Surgical Sets & Procedure Kits

- **Extraction Forceps**  
Elevators & Root Picks  
Periodontal Scalers & Curettes  
Restorative & Orthodontic Tools  
Diagnostic Sets & Mirrors  
Endodontic Instruments  
Impression Trays & Accessories

- **Bone Levers & Osteotomes**  
Bone Cutting Forceps & Rongeurs  
Chisels, Gouges & Mallets  
Plates, Screws & Fixation Tools  
Orthopedic Retractors & Clamps  
Screwdrivers, Taps & Drill Guides

- **Speculums & Retractors**  
Vaginal & Uterine Dilators  
Forceps & Scissors  
Obstetric Delivery Instruments  
Diagnostic & Examination Sets

- **Ear, Nose & Throat Diagnostic Sets**  
Nasal Speculums & Forceps  
Laryngoscopy & Tonsil Instruments  
Eye Speculums, Scissors & Forceps  
Micro Ophthalmic Instruments

- **Single-Use Scalpels & Forceps**  
Sterile Disposable Kits  
Plastic Handle Instruments  
Ready-to-Use Procedure Packs

- **General Veterinary Surgery Sets**  
Dental & Orthopedic Tools for Animals  
Diagnostic & Restraint Instruments

Precision. Partnership. Progress.

## Scalers, supragingival

Scaler, supragingivaler Einsatz

Limpiadores de sarro, supragingival

Instruments à détartrer, supragingival

Ablatori di tartaro, supragingival


## Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia



**SP29-116-01**

Fig. 1

**SP29-116-02**

Fig. 2

**SP29-116-03**

Fig. 3

**SP29-116-32**

Fig. 32

**SP29-116-33**

Fig. 33


**SP29-116-34**

Fig. 34


**SP29-116-01 -  
SP29-116-34**  
17,5 cm / 6 7/8"  
solid handle, 6 mm




**SP29-117-33 -  
 SP29-121-07**  
 17,5 cm / 6 7/8"  
 solid handle, 6 mm




**Jaquette  
 SP29-117-33**  
 Fig. 33




**Jaquette  
 SP29-117-34**  
 Fig. 34




**Jaquette  
 SP29-117-35**  
 Fig. 35




**Cumine  
 SP29-118-04**  
 Fig. 33




**Mitchell  
 SP29-119-04**  
 Fig. 34




**SP29-120-33**  
 Fig. 33/34



**Younger-Good  
 SP29-121-07**  
 Fig. 33/34




2,6 mm



1,6 mm



4,2 mm



1,0 mm

## Scalers, supragingival

Scaler, supragingivaler Einsatz  
 Limpiadores de sarro, supragingival  
 Instruments à détartrer, supragingival  
 Ablatori di tartaro, supragingival

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia



**Towner**  
**SP29-122-05**

Fig. 5/33



Incisors, premolars,  
 universal, also approximal



**McCall**  
**SP29-122-13**

Fig. 13S/14S



Lateral teeth, mesial/distal  
 Particularly useful for working in  
 interdental spaces in molar region




**McCall**  
**SP29-122-14**

Fig. 13SM/14



Lateral teeth; mesial/distal,  
 like 13S/14S, with rounded tip



**Taylor**  
**SP29-122-15**

Fig. T 2/3



Incisors, premolars, molars,  
 universal  
 Particularly useful for working  
 in interdental spaces in molar  
 region

**SP29-122-05 -**

**SP29-122-15**

17,5 cm / 6 7/8"

hollow handle, 8 mm

GR = grey

Y = yellow

R = red


G = green

V = violet

W = white

B = blue


N = black



**Towner**  
**SP29-122-16**  
 Fig. 15/30




For removing stubborn  
 tartar/scale



**SP29-122-17**  
 Fig. U 135




for hard-to-access  
 approximal surfaces



**SP29-122-18**  
 Fig. Cl 2/3




Incisors, premolars, molars,  
 universal



**SP29-122-19**  
 Fig. M 23



Incisors, premolars,  
 mesial/distal



**Mc Call**  
**SP29-123-18**  
 Fig. 17/18



universal


**SP29-122-16 -**  
**SP29-123-18**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm

## Scalers, supragingival

Scaler, supragingivaler Einsatz  
 Limpiadores de sarro, supragingival  
 Instruments à détartrer, supragingival  
 Ablatori di tartaro, supragingival

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia




**SP29-124-06 -  
 SP29-125-02**  
 17,5 cm / 6 3/4"  
 hollow handle, 8 mm

GR = grey    Y = yellow    R = red    G = green    V = violet    W = white    B = blue    N = black



1/4





**Zürich**  
**SP29-126-17**



For cleaning cavities  
 following apicœctomy



**SP29-127-17**  
 Fig. CC 6



**SP29-128-17**  
 Fig. T 2/3



Incisors, premolars,  
 molars, universal


**SP29-126-17-**  
**SP29-128-17**  
 17,5 cm / 6 3/4"  
 hollow handle, 8 mm

## Scalers, supragingival

Scaler, supragingivaler Einsatz  
Limpiadores de sarro, supragingival  
Instruments à détartrer, supragingival  
Ablatori di tartaro, supragingival



## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia




**Jaquette  
SP29-129-32**

Fig. 31/32



**Jaquette  
SP29-129-33**

Fig. 30/33



**Jaquette  
SP29-129-35**


Fig. 34/35




**SP29-129-32 -  
SP29-129-35**  
17,5 cm / 6 7/8"  
hollow handle, 8 mm




**SP29-130-31 -  
 SP29-130-36**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm





**Goldman-Fox  
 SP29-130-31**  
 Fig. 1





**Goldman-Fox  
 SP29-130-32**  
 Fig. 2





**Goldman-Fox  
 SP29-130-33**  
 Fig. 3






**Goldman-Fox  
 SP29-130-34**  
 Fig. 4



**Goldman-Fox  
 SP29-130-35**  
 Fig. 5



**Goldman-Fox  
 SP29-130-36**  
 Fig. 6



## Scalers, supragingival, bone files

Scaler, supragingivaler Einsatz, Knochenfeilen  
 Limpiadores de sarro, supragingival, limas de hueso  
 Instruments à détartrer, supragingival, râpes à os  
 Ablatori di tartaro, supragingival, lime e raspe per osso

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia



**Hirschfeld  
 SP29-131-03**

Fig. 3/7



**Hirschfeld  
 SP29-131-05**

Fig. 5/11




File

For removal of granulosomatous tissue from supra- and infrabone pockets


1/4

**SP29-131-03 -  
 SP29-131-05**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm



**Mercitan  
 SP29-132-05**

Fig. SH 5




**Ward-Kirkland  
 SP29-132-10**

Fig. 1



**Crane-Kaplan  
 SP29-132-20**

Fig. 2




**Crane-Kaplan  
 SP29-132-30**

Fig. 3

1/4


**SP29-132-05 -  
 SP29-132-30**  
 17,5 cm / 6 7/8"  
 solid handle, 6 mm




**White  
 SP29-133-12**

Fig. 12B


1/4





**SP29-134-01**




modified



**SP29-134-02**



**Orban**  
**SP29-134-06**  
 Fig. 12/13




**SP29-134-08**  
 Fig. 12/12



**SP29-134-01 -**  
**SP29-134-08**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm

For removal of granulosomatous tissue from supra- and infrabone pockets




**Periodontal file**  
Parodontalfeile  
Lima parodontal  
Lime parodontale  
Lima periodontale


**Scalers, periodontology**  
Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia




**Roncati**  
**SP29-135-00**  
17,5 cm / 6 7/8"  
hollow handle, 10 mm



**Roncati**  
**SP29-135-01**  
17,5 cm / 6 7/8"  
hollow handle, 10 mm



**Rhodes**  
**SP29-136-01**  
Fig. 36/37  
17,5 cm / 6 7/8"  
hollow handle, 8 mm



**Ochsenein  
 SP29-137-03**

Fig. 3

**Ochsenein  
 SP29-137-04**

Fig. 4

**Fedi  
 SP29-138-01**

Fig. 1

**Fedi  
 SP29-138-02**

Fig. 2


**Fedi  
 SP29-138-03**

Fig. 3

For flap operation. Also for minor bone corrections

**SP29-137-03 -  
 SP29-138-03**


17,5 cm / 6 7/8"  
 hollow handle, 8 mm



**Columbia**  
**SP29-139-13**  
 Fig. 13/14




Incisors and lateral teeth,  
 universal



**Columbia**  
**SP29-139-20**  
 Fig. 2R/2L




Incisors, universal




**Columbia**  
**SP29-139-40**  
 Fig. 4R/4L




Lateral teeth, universal



**SP29-140-02**  
 Fig. GX4



Incisors and lateral teeth,  
 universal



**Mc Call**  
**SP29-141-02**  
 Fig. 17S/18S



universal  
 bifurcate root surfaces

**SP29-139-13 -**  
**SP29-141-02**  
 17,5 cm / 6 7/8"  
 Universal curette  
 hollow handle, 8 mm


- GR grey
- Y yellow
- R red
- G green
- V violet
- W white
- B blue
- C light-blue

## Curettes, subgingival

Küretten, subgingivaler Einsatz  
Curetas, subgingival  
Curettes, subgingival  
Curette, subgingivale


## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia




GR

1/4



N

1/4



Y

1/4

	Langer	Langer	Langer
St	SP29-142-02	SP29-142-04	SP29-142-06
Ti	SP29-143-02	SP29-143-04	SP29-143-06
	Fig. 1/2	Fig. 3/4	Fig. 5/6



Mandible, lateral teeth,  
mesial / distal



Maxilla, lateral teeth,  
mesial / distal



Maxilla and mandible,  
incisors, universal

### Langer SP29-142-02 - SP29-143-06

17,5 cm / 6 7/8"

Universal curette  
hollow handle, 8 mm

GR = grey

Y = yellow

R = red

G = green

V = violet

W = white

B = blue

N = black



**Roncati Scaler**  
**SP29-144-00**  
 17,5 cm / 6 7/8"  
 Scaler  
 hollow handle, 10 mm



**Roncati**  
**SP29-144-01**  
 17,5 cm / 6 7/8"  
 Curette, micro, straight  
 hollow handle, 10 mm



**Roncati**  
**SP29-144-02**  
 17,5 cm / 6 7/8"  
 Curette, micro, curved  
 hollow handle, 10 mm



**Roncati**  
**SP29-145-00**  
 17,5 cm / 6 7/8"  
 Curette, micro, curved  
 hollow handle, 10 mm



Why scaling instruments should always be sharp:


1. Sharp instruments allow precision work with little application of force.
2. Blunt instruments polish the dental stone, bruise the soft tissue and require more application of force.

What you need for instrument sharpening:

1. Sharpening stones (rough, fine, and round for burr removal)
2. Sharpening oil (prevents excessive heat build-up; preserves the sharpening capacity of the sharpening stone)
3. Test stick (made of acrylic, for sharpness test)
4. Good light (sharp edges do not reflect light)
5. Clamping device (e.g. vise)

What to check after sharpening:

1. Roundness of back surfaces
2. Continuous sharpness of cutting edges
3. Perfect shape of working end (e.g. check against catalog illustration)
4. After several sharpening cycles, verify that the instrument is still strong enough for application. Replace instrument if necessary.



Color coding of scalers and Gracey curettes

1. All scalers are color-coded green – 3 points distributed across the diameter. At the same time, the color code marks the finer working end.
2. The Gracey Standard, Gracey Rigid and Gracey Mini curettes also have the color code on the side of the finer working end.  
Example: Gracey 1/2 = color code applied on side of working end no. 1.
3. For even better differentiation between Gracey Standard and Gracey Mini on the one hand and Gracey Rigid on the other hand, Gracey Rigid curettes are coded with six color points.
4. Color-coded Gracey Standard curettes are also available with a 10 mm handle.


left: color  
(lower figure,  
e.g. insert 3)



right: no color  
(higher figure,  
e.g. insert 4)



	∅		
Gracey Mini	8 mm	x	
Gracey Standard	8 mm	x	
Gracey Big	10 mm	x	
Gracey Rigid	8 mm		x



Fields of application for Gracey curettes


-  1-2,3-4, 5-6 Incisors
-  7-8, 9-10 Premolars and molars (buccal/lingual)
-  11-12, 15-16 Premolars and molars, mesial  
Furcation: mesial portion of the distal root
-  13-14, 17-18 Premolars and molars, distal  
Furcation: distal portion of the mesial root

**Gracey mini curettes, subgingival**

Gracey Mini-Küretten, subgingivaler Einsatz  
 Gracey mini curetas, subgingival  
 Gracey mini curettes, subgingival  
 Gracey mini curette, subgingivale


**Scalers, periodontology**

Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia




**Gracey Mini**  
**SP29-146-02 -**  
**SP29-146-18**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm


GR = grey    Y = yellow    R = red    G = green    V = violet    W = white    B = blue    N = black




**SP29-146-10**  
 Fig. 9/10




Premolars, molars,  
 lingual/buccal,  
 first shaft longer,  
 working end shorter




**SP29-146-12**  
 Fig. 11/12




Premolars, molars,  
 mesial, mesiolingual,  
 mesiobuccal,  
 first shaft longer,  
 working end shorter




**SP29-146-14**  
 Fig. 13/14




Premolars, molars, distal,  
 first shaft longer,  
 working end shorter




**SP29-146-16**  
 Fig. 15/16



Premolars, molars, mesial,  
 first shaft longer,  
 working end shorter



**SP29-146-18**  
 Fig. 17/18




Premolars, molars, distal,  
 first shaft longer,  
 working end shorter

**Gracey standard curettes, subgingival**


Gracey Standard-Küretten, subgingivaler Einsatz  
 Gracey curetas estándar, subgingival  
 Gracey curettes standard, subgingival  
 Gracey curette standard, subgingivale

**Scalers, periodontology**


Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia



**SP29-147-02**  
Fig. 1/2




Incisors, universal



**SP29-147-04**  
Fig. 3/4




Incisors, premolars,  
universal



**SP29-147-06**  
Fig. 5/6



Incisors, premolars,  
universal




**SP29-147-08**  
Fig. 7/8




Premolars, molars,  
lingual / buccal

**Gracey Standard**  
**SP29-147-02 -**  
**SP29-147-18**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm






**SP29-147-10**  
Fig. 9/10




Premolars, molars,  
lingual / buccal




**SP29-147-12**  
Fig. 11/12




Premolars, molars, mesial,  
mesiolingual, mesiobuccal



**SP29-147-14**  
Fig. 13/14




Premolars, molars, distal




**SP29-147-16**  
Fig. 15/16



Premolars, molars, mesial



**SP29-147-18**  
Fig. 17/18




Premolars, molars, distal

**Gracey Rigid curettes, subgingival**


Gracey Rigid-Küretten, subgingivaler Einsatz  
 Gracey Rigid curetas, subgingival  
 Gracey Rigid curettes, subgingival  
 Gracey Rigid curette, subgingivale

**Scalers, periodontology**


Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia



**SP29-148-02**  
Fig. 1/2




Incisors, universal,  
extra rigid, for stubborn  
dental plaque



**SP29-148-04**  
Fig. 3/4




Incisors, premolars,  
universal,  
extra rigid, for stubborn  
dental plaque



**SP29-148-06**  
Fig. 5/6



Incisors, premolars,  
universal,  
extra rigid, for stubborn  
dental plaque



**SP29-148-08**  
Fig. 7/8




Premolars, molars, lingual,  
buccal,  
extra rigid, for stubborn  
dental plaque


**Gracey Rigid**  
**SP29-148-02 -**  
**SP29-148-18**

17,5 cm / 6 7/8"  
hollow handle, 8 mm


- GR grey
- Y yellow
- R red
- G green
- V violet
- W white
- B blue
- N black



**SP29-148-10**  
Fig. 9/10




Premolars, molars, lingual, buccal, extra rigid, for stubborn dental plaque




**SP29-148-12**  
Fig. 11/12




Premolars, molars, mesial, mesiolingual, mesiobuccal, extra rigid, for stubborn dental plaque



**SP29-148-14**  
Fig. 13/14




Premolars, molars, distal, extra rigid, for stubborn dental plaque




**SP29-148-16**  
Fig. 15/16



Premolars, molars, mesial, extra rigid, for stubborn dental plaque



**SP29-148-18**  
Fig. 17/18




Premolars, molars, distal, extra rigid, for stubborn dental plaque

**Gracey curettes with 10 mm handle**


Gracey Küretten mit 10 mm-Griff  
 Gracey curetas para mangos 10 mm  
 Gracey curettes avec manches à 10 mm  
 Gracey curette con manici 10 mm

**Scalers, periodontology**


Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia



**SP29-149-02**  
Fig. 1/2




Incisors, universal



**SP29-149-04**  
Fig. 3/4




Incisors, premolars, universal



**SP29-149-06**  
Fig. 5/6



Incisors, premolars, universal




**SP29-149-08**  
Fig. 7/8



Premolars, molars, lingual, buccal


**Big Gracey**  
**SP29-149-02 -**  
**SP29-149-18**  
 17,5 cm / 6 7/8"  
 hollow handle, 10 mm

GR grey    Y yellow    R red    G green    V violet    W white    B blue    N black




**SP29-149-10**  
Fig. 9/10

Premolars, molars, lingual, buccal




**SP29-149-12**  
Fig. 11/12

Premolars, molars, mesial, mesiolingual, mesiobuccal,




**SP29-149-14**  
Fig. 13/14

Premolars, molars, distal





**SP29-149-16**  
Fig. 15/16

Premolars, molars, mesial



**SP29-149-18**  
Fig. 17/18

Premolars, molars, distal



## Curettes, subgingival

Küretten, subgingivaler Einsatz  
Curetas, subgingival  
Curettes, subgingival  
Curette, subgingivale

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



**SP29-150-01**

Fig. 1




Furcation curette  
0,9 mm  
buccal / lingual



**SP29-150-02**

Fig. 2



Furcation curette  
3 mm  
buccal / lingual



**SP29-150-03**

Fig. 3



Furcation curette  
0,9 mm  
mesial / distal



**SP29-150-04**

Fig. 4



Furcation curette  
1,3 mm  
mesial / distal

**SP29-150-01 -  
SP29-150-04**

17,5 cm / 6 7/8"  
hollow handle, 8 mm




**SP29-150-05**

Fig. 5



Furcation curette  
 0,9 mm  
 labial / lingual



**Prichard**  
**SP29-151-01**

Fig. 1/2



Bone curette,  
 extra long, robust  
 construction, double  
 cutting

### Gracey mini curettes, subgingival

Gracey Mini-Küretten, subgingivaler Einsatz

Gracey mini curetas, subgingival

Gracey mini curettes, subgingival

Gracey mini curette, subgingivale

### Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia



**SP29-152-06**

Ti

Fig. 5/6

17,5 cm / 6 7/8"


hollow handle, 10 mm

0,8 mm

1/4

Incisors, premolars,  
universal,  
first shaft longer,  
working end shorter

Incisivi, premolari,  
universali,  
1° stelo più lungo,  
estremità di lavoro più breve



**SP29-152-08**

Ti

Fig. 7/8



17,5 cm / 6 7/8"

hollow handle, 10 mm

0,8 mm

1/4

Premolars, molars, lingual/buccal,  
first shaft longer,  
working end shorter



Premolars, molars,  
 mesial, mesiolingual,  
 mesiobuccal,  
 first shaft longer, working  
 end shorter

Premolari, molari,  
 mesiali, mesiolinguali,  
 mesiobuccali,  
 1° stelo più lungo, estremità  
 di lavoro più breve

1/4



**SP29-152-12**

Ti

Fig. 11/12

17,5 cm / 6 7/8"

hollow handle, 10 mm



Premolars, molars, distal,  
 first shaft longer, working  
 end shorter

1/4

**SP29-152-14**

Ti

Fig. 13/14

17,5 cm / 6 7/8"

hollow handle, 10 mm

### Gracey standard curettes, subgingival

Gracey Standard-Küretten, subgingivaler Einsatz

Gracey curetas estándar, subgingival

Gracey curettes standard, subgingival

Gracey curette standard, subgingivale


### Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia



**SP29-153-06**

Ti

Fig. 5/6

17,5 cm / 6 7/8"

hollow handle, 10 mm




**SP29-153-08**

Ti

Fig. 7/8

17,5 cm / 6 7/8"


hollow handle, 10 mm



Premolars, molars, mesial,  
 mesiolingual, mesiobuccal

**SP29-153-12**  
 Fig. 11/12  
 17,5 cm / 6 7/8"  
 hollow handle, 10 mm

Ti



Premolars, molars, distal

**SP29-153-14**  
 Fig. 13/14  
 17,5 cm / 6 7/8"  
 hollow handle, 10 mm

Ti

## Sharpening stones

Schleifsteine  
Piedras de afilar  
Meules  
Pietre per affilare

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



### Record Arkansas SP29-154-00

115 x 40 x 10/2 mm  
multi shape, white

For pregrinding of blunt instruments. The rounded edges and bevelled surfaces enable the sharpening of radius and surfaces



### Original Arkansas SP29-154-01

100 x 25 x 7 mm  
flat, grey-white

For repeated grinding and finishing of scalers and cures



### Original Arkansas SP29-154-02

100 mm lang, Ø 8 mm  
round, grey-white

For burr removing from cutting edges



### India-Korund SP29-154-05

115 x 45 x 2 - 10 mm

delicate, brown-beige  
Para preafilar y volverles a dar forma a instrumentos



### India-Korund SP29-154-06

100 mm lang, Ø 5 - 10 mm

round, delicate, brown-beige  
For burr removing from cutting edges



### SP29-154-08

75 mm lang, Ø 6,0 mm

Test stick plastic  
package of 6 pieces  
For testing the sharpness of scalers and cures



**Record Arkansas**  
**SP29-155-00**  
150 x 50 mm  
in wooden case




**SP29-154-09**  
Sharpening stone oil

## Periodontal probes

Parodontometer  
Sondas parodontales  
Sondes parodontales  
Sonde parodontali

## Scalers, periodontology


Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



**SP29-156-01**



**SP29-157-40**



**SP29-158-08**


Fig. 8

Grad. 3-6-8-11

**SP29-158-10**

Fig. 10

Grad. 1-2-3-5-7-8-9-10




### Iglhaut

**SP29-156-01**

17,5 cm / 6 7/8"

solid handle, 6 mm




### Nabers

**SP29-157-40**

17,5 cm / 6 7/8"

hollow handle, 8 mm

Furcation probe, graduated working end



**SP29-158-11**

Fig. 11

Grad. 3-6-8-11

**SP29-159-00**


Fig. W

Grad. 1-2-3-5-7-8-9-10

### SP29-158-08 - SP29-159-00

16 cm / 6 1/8"

solid handle, 6 mm



<b>ST</b>	<b>SP29-160-11</b>	<b>SP29-160-13</b>	<b>SP29-160-21</b>	<b>SP29-160-23</b>
<b>TI</b>	<b>SP29-161-11</b>	<b>SP29-161-13</b>	<b>SP29-161-21</b>	<b>SP29-161-23</b>
	Fig. CPG 11,5 WHO	Fig. CPG 12	Fig. CP 15 UNC	Fig. CPNG 22
	Grad. 3,5-5,5-8,5-11,5 with ball	Grad. 3-6-9-12	Grad. 1-2-3-4-...-14-15	Grad. 2-4-6-8-10-12




Probe handles see page 244

**Periodontal probes**

Parodontometer  
 Sondas parodontales  
 Sondes parodontales  
 Sonde parodontali


**Scalers, periodontology**

Reinigungsinstrumente, Parodontologie  
 Limpiadores de sarro, parodontología  
 Instruments à détartrer, parodontologie  
 Ablatori di tartaro, parodontologia




**Merrit A**  
**SP29-162-00**

Grad. 1-2-3-5-7-8-9-10




**Merrit E**  
**SP29-162-01**

Grad. 1-2-3-5-7-8-9-10




**Goldman-Fox**  
**SP29-162-04**

Grad. 1-2-3-5-7-8-9-10




**Michigan**  
**SP29-162-06**

Grad. 3-6-8




**Mercury**  
**SP29-162-08**

Grad. 3-6-9-12




**SP29-162-10**

Fig. CPNG 12S  
 Grad. 3-6-9-12




**WHO**  
**SP29-162-11**

Fig. CPG 11,5 WHO  
 Grad. 3,5-5,5-8,5-11,5




**SP29-162-12**

Fig. CPG 12  
 Grad. 3-6-9-12




**SP29-162-15**

Fig. CP 15 UNC  
 Grad. 1-2-3-4-5 ... 14-15



**SP29-162-22**

Fig. CPNG 22  
 Grad. 2-4-6-8-10-12



**SP29-162-26**
















Fig. 26G  
 Grad. 1-2-3-5-7-8-9-10

**SP29-162-00 -**

**SP29-162-26**

17,5 cm / 6 7/8"  
 solid handle, 6 mm



	 <p>0,8 mm</p> <p>1/4</p>	 <p>0,6 mm</p>	 <p>0,6 mm</p>	 <p>0,5 mm</p> <p>0,8 mm</p>	 <p>0,5 mm</p> <p>0,8 mm</p>
	<p><b>Williams-Goldman-Fox</b>  <b>SP29-163-02</b>                  Grad. 1-2-3-5-7-8-9-10</p>	<p><b>SP29-163-04</b>                  Fig. CP 11/23                  Grad. 3-6-8-11</p>	<p><b>SP29-163-06</b>                  Fig. CP 12/23                  Grad. 3-6-9-12</p>	<p><b>SP29-163-08</b>                  Fig. GF/3                  Grad. 1-2-3-5-7-8-9-10</p>	<p><b>SP29-163-10</b>                  Fig. UM/3                  Grad. 3-6-8</p>
	 <p>0,5 mm</p> <p>0,8 mm</p> <p>1/4</p>	 <p>0,5 mm</p>	 <p>0,5 mm</p>	 <p>0,9 mm</p>	 <p>0,6 mm</p>
	<p><b>SP29-163-12</b>                  Fig. W/3                  Grad. 1-2-3-5-7-8-9-10</p>	<p><b>SP29-163-14</b>                  Fig. GF/6                  Grad. 1-2-3-5-7-8-9-10</p>	<p><b>SP29-163-16</b>                  Fig. UM/6                  Grad. 3-6-8</p>	<p><b>SP29-163-18</b>                  Fig. W/6                  Grad. 1-2-3-5-7-8-9-10</p>	<p><b>SP29-163-20</b>                  Fig. GF/17                  Grad. 1-2-3-5-7-8-9-10</p>
<p><b>SP29-163-02 -</b>  <b>SP29-163-20</b>                  17,5 cm / 6 7/8"                  hollow handle, 8 mm</p>	 <p>0,5 mm</p>	 <p>0,9 mm</p>	 <p>0,6 mm</p>	 <p>0,6 mm</p>	 <p>0,9 mm</p>

## Periodontal probes

Parodontometer  
Sondas parodontales  
Sondes parodontales  
Sonde parodontali

## Scalers, periodontology

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



0,5 mm




1/4

### SP29-163-22

Fig. UM/17

Grad. 3-6-8



0,5 mm

0,5 mm



### SP29-163-24

Fig. W/17

Grad. 1-2-3-5-7-8-9-10



0,5 mm


0,9 mm



### SP29-163-26

Fig. GF/23

Grad. 1-2-3-5-7-8-9-10



0,6 mm


0,6 mm



### SP29-163-28


Fig. UM/23

Grad. 3-6-8



0,5 mm


0,6 mm



### SP29-163-30

Fig. W/23

Grad. 1-2-3-5-7-8-9-10



0,5 mm


### SP29-163-22 - SP29-163-30

17,5 cm / 6 7/8"

hollow handle, 8 mm




**Cross**  
**SP29-164-10**  
17,5 cm / 6 7/8"  
hollow handle, 8 mm



**Zürich**  
**SP29-165-02**  
17,5 cm / 6 7/8"  
hollow handle, 8 mm

-  GR grey
-  Y yellow
-  R red
-  G green
-  V violet
-  W white
-  B blue
-  C light-blue


**FinoLine**



**SP29-165-03**  
Fig. 3




**SP29-165-06**  
Fig. 6



**SP29-165-08**  
Fig. 8 A




**SP29-165-09**  
Fig. 9



**SP29-165-12**  
Fig. 12




**SP29-165-17**  
Fig. 17




**SP29-165-23**  
Fig. 23



**SP29-165-30**  
Fig. 3A




**SP29-165-60**  
Fig. 6XL



**White**  
**SP29-165-12w**  
Fig. 12


**SP29-165-03 -**  
**SP29-165-12w**  
17,5 cm / 6 7/8"  
extra delicate tooth probe  
solid handle, 6 mm



**SP29-166-03**

Fig. 3


Detection of furcations and for easy entry into interproximal spaces



**SP29-166-06**

Fig. 6

Detection of plaque, caries, root canals and furcations



**SP29-166-08**

Fig. 8A




**SP29-166-09**

Fig. 9



**Briault**  
**SP29-166-11**

Fig. EXD 11




**SP29-166-12**

Fig. EXD 11-12




**SP29-166-03 -**  
**SP29-166-12**  
17,5 cm / 6 7/8"  
solid handle, 6 mm




**White**  
**SP29-166-17**  
 Fig. 17

Mainly for sub-gingival examination and approx. caries




**White**  
**SP29-166-23**  
 Fig. 23


For general examination, especially for incisors



**SP29-166-33**  
 Fig. 3A



**SP29-166-54**  
 Fig. 54




**SP29-166-66**  
 Fig. 6XL

Detection of plaque, caries, root canals and furcations

**SP29-166-17 -**  
**SP29-166-66**  
 17,5 cm / 6 7/8"  
 solid handle, 6 mm




**SP29-167-01 -  
SP29-167-23**  
17,5 cm / 6 7/8"  
hollow handle, 8 mm




**SP29-167-01**  
Fig. UT 1/2




**SP29-167-02**  
Fig. UT 2/2




**SP29-167-03**  
Fig. 3CH/3CH




**SP29-167-04**  
Fig. 6/3




**SP29-167-06**  
Fig. 23/6







**SP29-167-16**  
Fig. 16/17




**SP29-167-17**  
Fig. 17/23




**SP29-167-23**  
Fig. 23/17A






**Nabers**  
**SP29-168-01**  
 Fig. 1







**Nabers**  
**SP29-168-02**  
 Fig. 2




**SP29-168-05**  
 Fig. 5/6



**SP29-168-23**  
 Fig. 16/23



**SP29-168-01 -**  
**SP29-168-23**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm



**SP29-169-08**  
 Fig. EXG 8


**SP29-169-09**  
 Fig. EXG 9

**SP29-169-10**  
 Fig. EXG 10


**SP29-169-23**  
 Fig. EXG 23



Probe handles




**SP29-170-01**  
 Fig. GM 1A




**SP29-170-02**  
 Fig. GM 2A



**SP29-171-00**  
 Fig. GM 0




**SP29-171-01**  
 Fig. GM 1



**SP29-170-01 -  
 SP29-171-34**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm  
 double ended

- GR  
 grey
- Y  
 yellow
- R  
 red
- G  
 green
- V  
 violet
- W  
 white
- B  
 blue
- C  
 light-blue




**SP29-171-02**  
 Fig. GM 2



**SP29-171-03**  
 Fig. GM 3




**SP29-171-04**  
 Fig. GM 4




**Zürich**  
**SP29-171-34**  
 Ø 1,9 mm







**SP29-172-00 -  
 SP29-172-44**  
 17,5 cm / 6 7/8"  
 hollow handle, 8 mm  
 double ended




**SP29-172-24**  
 Fig. GM 17



**SP29-172-26**  
 Fig. GM 18



**SP29-172-28**  
 Fig. GM 19



**SP29-172-30**  
 Fig. GM 17W



**SP29-172-32**  
 Fig. GM 18W




**SP29-172-34**  
 Fig. GM 19W




**SP29-172-35**  
 Fig. GM 61/62




**SP29-172-36**  
 Fig. GM 63/64




**SP29-172-38**  
 Fig. GM 65/66






**SP29-172-40**  
 Fig. GM 31L



**SP29-172-42**  
 Fig. GM 32L



**SP29-172-44**  
 Fig. GM 33L




**Enamel chisels**

Schmelzmeißel  
Cinceles para esmalte  
Ciseaux à émail  
Tagliasmalto

**Scalers, periodontology**

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



**SP29-173-01**  
Fig. 1




**SP29-173-02**  
Fig. 2




**SP29-173-26**  
Fig. 26



**SP29-173-27**  
Fig. 27




**SP29-173-28**  
Fig. 28



**SP29-173-29**  
Fig. 29




**SP29-173-52**  
Fig. 51/52




**SP29-173-54**  
Fig. 53/54




**SP29-173-64**  
Fig. 63/64



**SP29-173-66**  
Fig. 65/66



**SP29-173-78**  
Fig. 77/78



**SP29-173-80**  
Fig. 79/80



1/4


**Black**

**SP29-173-01 -  
SP29-173-80**


17,5 cm / 6 7/8"  
solid handle, 6 mm  
double ended




1/4






**SP29-174-02**  
Fig. 1/2



**SP29-174-04**  
Fig. 3/4



**SP29-174-06**  
Fig. 5/6



solid handle, 6 mm  
double ended

**Wedelstaedt**  
**SP29-174-02 -**  
**SP29-174-06**  
17,5 cm / 6 7/8"

**Root canal instruments**

Nervkanal-Instrumente  
Instrumentos para canales  
Instruments à canaux  
Strumenti per trattamenti canalari

**Scalers, periodontology**

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



1/4


**SP29-175-31 -**

**SP29-176-01**

17,5 cm / 6 7/8"


hollow handle, 8 mm

double ended




**Zürich**  
**SP29-175-31**

Fig. 30/40




**Zürich**  
**SP29-175-51**


Fig. 50/90



**Stewart**  
**SP29-176-01**



Plugger



Plugger



Endo probe



grey



yellow



red



green



violet



white



blue




black




1/4


**SP29-177-03 -  
 SP29-177-12**  
 17,5 cm / 6 7/8"  
 solid handle, 6 mm  
 single ended




**SP29-177-03**  
 Fig. 30




**SP29-177-04**  
 Fig. 40




**SP29-177-05**  
 Fig. 50




**SP29-177-06**  
 Fig. 60



**SP29-177-09**  
 Fig. 90



**SP29-177-10**  
 Fig. 100



**SP29-177-12**  
 Fig. 120


Plugger

**Root canal instruments**

Nervkanal-Instrumente  
Instrumentos para canales  
Instruments à canaux  
Strumenti per trattamenti canalari


**Scalers, periodontology**

Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia




**SP29-178-01 -  
SP29-179-06**  
16 cm / 6 1/8"  
solid handle, 6 mm  
single ended


- GR grey
- Y yellow
- R red
- G green
- V violet
- W white
- B blue
- N black



**SP29-180-37**  
 Fig. 7




**SP29-180-38**  
 Fig. 8




**SP29-180-39**  
 Fig. 9



**SP29-180-40**  
 Fig. 10



**SP29-180-41**  
 Fig. 11




**SP29-180-37 -  
 SP29-180-41**  
 17,5 cm / 6 7/8"  
 solid handle, 6 mm  
 double ended

## Gingivectomy knives



Gingivektomiemesser  
Bisturías para la gingivectomía  
Bistouris pour la gingivectomie  
Bisturie per gengivectomia

## Scalers, periodontology


Reinigungsinstrumente, Parodontologie  
Limpiadores de sarro, parodontología  
Instruments à détartrer, parodontologie  
Ablatori di tartaro, parodontologia



SP29-181-03



SP29-181-04




**Allen-Iglhaut**

**SP29-181-03 -**

**SP29-181-04**


16 cm / 6 1/8"

hollow handle, 10 mm  
double ended




**Kirkland**  
**SP29-182-51**


For gingivectomy or gingivoplastic surge  
 Very effective in retro molar region



**Kirkland**  
**SP29-182-61**




**Kirkland**  
**SP29-183-61**




**Orban**  
**SP29-184-01**

Fig. 1



**Orban**  
**SP29-184-02**


Fig. 2



**Orban**  
**SP29-184-12**

Fig. 1/2

Counter angled shank, especially suitable for excision of  
 interproximal tissue in gingivectomy



**SP29-182-51 -**  
**SP29-184-12**  
 16 cm / 6 1/8"  
 solid handle, 6 mm  
 single ended

# Integrity in Every Incision.

Thank you for taking the time to learn about Surgiplex SRL – your trusted partner in surgical innovation and healthcare supply.

**Head Office**

Strada Vanatori 3, Pravalia 6,  
Sector 5, Bucharest, Romania

**Email:** [surgiplexsrl@gmail.com](mailto:surgiplexsrl@gmail.com)

**WhatsApp:** +92 333 8635932

**EU Contact:** +40 765 651 717

**GCC Contact:** +974 3397 9907

**Website:** [www.surgiplexsrl.com](http://www.surgiplexsrl.com)

# SURGIPILEX