

SPX-29

SCALERS, PERIODONTOLOGY

Efficiency. Precision
Care — gentle removal,
superior control.

Surgiplex. The Difference You Hold



www.surgiplexsrl.com

ABOUT US

The Bridge Between Manufacturing Excellence and Global Supply

Surgiplex SRL is a Romania-based export and distribution company specializing in the global supply of high-quality surgical instruments. Operating from Bucharest, Surgiplex acts as the commercial and logistics hub of the Surgilabs Group, a partnership designed to unify manufacturing expertise with European-level quality assurance and distribution.

In collaboration with Surgilabs (Pakistan) — our dedicated manufacturing partner — we ensure a seamless connection between production precision and international healthcare delivery. Surgilabs is responsible for crafting each surgical instrument to global standards, while Surgiplex manages export, compliance, and customer relations, ensuring every client receives products that reflect reliability, safety, and performance.

OUR COLLABORATION TRANSFORMS PRECISION INTO RELIABILITY — EVERY INSTRUMENT WE EXPORT CARRIES A PROMISE OF PERFORMANCE.

Together, Surgiplex and Surgilabs form a vertically integrated network — one that transforms local craftsmanship into a global benchmark of quality and care.



VISION | MISSION | CORE VALUES

To be a trusted global partner in delivering innovative, precision-crafted surgical instruments that improve healthcare outcomes worldwide.

To connect healthcare providers with superior-quality surgical products through ethical practices, seamless logistics, and exceptional customer care, powered by our dedicated production partner, Surgilabs.

- INTEGRITY

Transparent and ethical in every process.

- EXCELLENCE

Instruments that exceed international standards.

- INNOVATION

Modern designs and flexible solutions.

- ACCESSIBILITY

Affordable instruments for healthcare worldwide.

- PARTNERSHIP

Building long-term networks with distributors and providers.

20+

Active distributors and partners worldwide.

5-7%

Annual revenue investment on R&D

Precision. Partnership. Progress.

PRODUCT RANGE

From Design to Delivery — A Complete Range of Surgical Solutions

At Surgiplex SRL, our product range reflects a commitment to quality, innovation, and precision across every discipline of surgery.

In partnership with Surgilabs (Pakistan), we manufacture and supply a comprehensive line of surgical, dental, and specialty instruments, engineered for performance, reliability, and ergonomics.

■ Precision:

Every edge, hinge, and grip is engineered to meet surgical accuracy.

■ Performance:

Designed for longevity, balance, and ease of handling.

■ Safety:

Manufactured under strict hygiene and sterilization protocols.

■ Innovation:

Continuous design refinement based on global medical feedback.



GENERAL SURGERY INSTRUMENTS
DENTAL INSTRUMENTS
ORTHOPEDIC INSTRUMENTS
GYNECOLOGY & OBSTETRICS
ENT & OPHTHALMIC INSTRUMENTS
SINGLE-USE INSTRUMENTS
CUSTOM & OEM INSTRUMENTS

100+ Total Categories

- **General Surgery Instruments**
Dissecting & Dressing Forceps
Surgical Scissors & Shears
Needle Holders & Clamps
Retractors (Self-Retaining & Hand-Held)
Hemostatic & Artery Forceps
Towel Clips & Sponge Forceps
Operating & Tissue Forceps
Surgical Sets & Procedure Kits

- **Extraction Forceps**
Elevators & Root Picks
Periodontal Scalers & Curettes
Restorative & Orthodontic Tools
Diagnostic Sets & Mirrors
Endodontic Instruments
Impression Trays & Accessories

- **Bone Levers & Osteotomes**
Bone Cutting Forceps & Rongeurs
Chisels, Gouges & Mallets
Plates, Screws & Fixation Tools
Orthopedic Retractors & Clamps
Screwdrivers, Taps & Drill Guides

- **Speculums & Retractors**
Vaginal & Uterine Dilators
Forceps & Scissors
Obstetric Delivery Instruments
Diagnostic & Examination Sets

- **Ear, Nose & Throat Diagnostic Sets**
Nasal Speculums & Forceps
Laryngoscopy & Tonsil Instruments
Eye Speculums, Scissors & Forceps
Micro Ophthalmic Instruments

- **Single-Use Scalpels & Forceps**
Sterile Disposable Kits
Plastic Handle Instruments
Ready-to-Use Procedure Packs

- **General Veterinary Surgery Sets**
Dental & Orthopedic Tools for Animals
Diagnostic & Restraint Instruments

Precision. Partnership. Progress.

Scalers, supragingival

Scaler, supragingivaler Einsatz

Limpiadores de sarro, supragingival

Instruments à détartrer, supragingival

Ablatori di tartaro, supragingival


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología


Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-116-01

Fig. 1




SP29-116-02

Fig. 2




SP29-116-03

Fig. 3




SP29-116-32

Fig. 32



SP29-116-33

Fig. 33



SP29-116-34

Fig. 34

SP29-116-01 -


SP29-116-34

17,5 cm / 6 7/8"


solid handle, 6 mm




**SP29-117-33 -
 SP29-121-07**
 17,5 cm / 6 7/8"
 solid handle, 6 mm




**Jaquette
 SP29-117-33**
 Fig. 33




**Jaquette
 SP29-117-34**
 Fig. 34




**Jaquette
 SP29-117-35**
 Fig. 35




**Cumine
 SP29-118-04**
 Fig. 33





**Mitchell
 SP29-119-04**
 Fig. 34



SP29-120-33
 Fig. 33/34



**Younger-Good
 SP29-121-07**
 Fig. 33/34



Scalers, supragingival

Scaler, supragingivaler Einsatz

Limpiadores de sarro, supragingival

Instruments à détartrer, supragingival

Ablatori di tartaro, supragingival

Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia



Towner
SP29-122-05

Fig. 5/33



Incisors, premolars,
universal, also approximal



McCall
SP29-122-13

Fig. 13S/14S



Lateral teeth, mesial/distal
Particularly useful for working in
interdental spaces in molar region




McCall
SP29-122-14

Fig. 13SM/14



Lateral teeth; mesial/distal,
like 13S/14S, with rounded tip



Taylor
SP29-122-15

Fig. T 2/3



Incisors, premolars, molars,
universal
Particularly useful for working
in interdental spaces in molar
region

SP29-122-05 -
SP29-122-15

17,5 cm / 6 7/8"

hollow handle, 8 mm

GR = grey

Y = yellow

R = red


G = green

V = violet

W = white

B = blue


N = black



Towner
SP29-122-16
Fig. 15/30



For removing stubborn
tartar/scale



SP29-122-17
Fig. U 135



for hard-to-access
approximal surfaces



SP29-122-18
Fig. CI 2/3




Incisors, premolars, molars,
universal



SP29-122-19
Fig. M 23



Incisors, premolars,
mesial/distal



Mc Call
SP29-123-18
Fig. 17/18



universal

SP29-122-16 -
SP29-123-18

17,5 cm / 6 7/8"

hollow handle, 8 mm

Scalers, supragingival

Scaler, supragingivaler Einsatz

Limpiadores de sarro, supragingival

Instruments à détartrer, supragingival

Ablatori di tartaro, supragingival


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología


Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




Hygienist
SP29-124-06

Fig. H6/H7



Incisivos, premolares,
mesial/distal
Esp. indicado para los
espacios interdentes



Hygienist
SP29-124-10

Fig. H6/H7



Incisivos, premolares,
mesial/distal
Esp. indicado para los
espacios interdentes



SP29-125-00

Fig. 204



Molars, universal



SP29-125-01

Fig. 204 S



Molars, universal



SP29-125-02

Fig. 204 SD



Incisors, premolars,
mesial/distal

SP29-124-06 -

SP29-125-02

17,5 cm / 6 3/4"

hollow handle, 8 mm

GR = grey

Y = yellow

R = red


G = green

V = violet

W = white

B = blue


N = black



Zürich
SP29-126-17



For cleaning cavities
following apicœctomy



SP29-127-17
Fig. CC 6



SP29-128-17
Fig. T 2/3



Incisors, premolars,
molars, universal

**SP29-126-17-
SP29-128-17**
17,5 cm / 6 3/4"
hollow handle, 8 mm

Scalers, supragingival

Scaler, supragingivaler Einsatz

Limpiadores de sarro, supragingival

Instruments à détartrer, supragingival

Ablatori di tartaro, supragingival


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología


Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




**Jaquette
SP29-129-32**

Fig. 31/32



**Jaquette
SP29-129-33**

Fig. 30/33



**Jaquette
SP29-129-35**

Fig. 34/35



SP29-129-32 -


SP29-129-35

17,5 cm / 6 7/8"

hollow handle, 8 mm



**SP29-130-31 -
SP29-130-36**
17,5 cm / 6 7/8"
hollow handle, 8 mm




**Goldman-Fox
SP29-130-31**

Fig. 1




**Goldman-Fox
SP29-130-32**

Fig. 2




**Goldman-Fox
SP29-130-33**

Fig. 3




**Goldman-Fox
SP29-130-34**

Fig. 4



**Goldman-Fox
SP29-130-35**

Fig. 5



**Goldman-Fox
SP29-130-36**

Fig. 6



Scalers, supragingival, bone files

Scaler, supragingivaler Einsatz, Knochenfeilen
 Limpiadores de sarro, supragingival, limas de hueso
 Instruments à détartrer, supragingival, râpes à os
 Ablatori di tartaro, supragingival, lime e raspe per osso

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
 Limpiadores de sarro, parodontología
 Instruments à détartrer, parodontologie
 Ablatori di tartaro, parodontologia



**Hirschfeld
SP29-131-03**

Fig. 3/7




**Hirschfeld
SP29-131-05**

Fig. 5/11




File

For removal of granulosomatous tissue from
supra- and infrabone pockets



**Mercitan
SP29-132-05**

Fig. SH 5



**Ward-Kirkland
SP29-132-10**

Fig. 1



**Crane-Kaplan
SP29-132-20**

Fig. 2



**Crane-Kaplan
SP29-132-30**

Fig. 3




**White
SP29-133-12**


Fig. 12B

**SP29-131-03 -
SP29-131-05**
17,5 cm / 6 7/8"
hollow handle, 8 mm


**SP29-132-05 -
SP29-132-30**
17,5 cm / 6 7/8"
solid handle, 6 mm





SP29-134-01



modified




SP29-134-02



Orban
SP29-134-06

Fig. 12/13



SP29-134-08

Fig. 12/12




SP29-134-01 -

SP29-134-08

17,5 cm / 6 7/8"


hollow handle, 8 mm

For removal of granulosomatous tissue from supra- and infrabone pockets



Periodontal file
 Parodontalfeile
 Lima parodontal
 Lime parodontale
 Lima periodontale


Scalers, periodontology
 Reinigungsinstrumente, Parodontologie
 Limpiadores de sarro, parodontología
 Instruments à détartrer, parodontologie
 Ablatori di tartaro, parodontologia




Roncati
SP29-135-00
 17,5 cm / 6 7/8"
 hollow handle, 10 mm



Roncati
SP29-135-01
 17,5 cm / 6 7/8"
 hollow handle, 10 mm



Mini chisel




Rhodes
SP29-136-01

Fig. 36/37


17,5 cm / 6 7/8"

hollow handle, 8 mm



For flap operation. Also for minor bone corrections

**SP29-137-03 -
 SP29-138-03**
 17,5 cm / 6 7/8"
 hollow handle, 8 mm




Columbia
SP29-139-13

Fig. 13/14



Incisors and lateral teeth,
universal




Columbia
SP29-139-20

Fig. 2R/2L



Incisors, universal



Columbia
SP29-139-40

Fig. 4R/4L




Lateral teeth, universal




SP29-140-02

Fig. GX4



Incisors and lateral teeth,
universal



Mc Call
SP29-141-02

Fig. 17S/18S



universal
bifurcate root surfaces

**SP29-139-13 -
SP29-141-02**

17,5 cm / 6 7/8"

Universal curette

hollow handle, 8 mm



GR grey



Y yellow



R red



G green



V violet



W white



B blue



C light-blue

Curettes, subgingival

Küretten, subgingivaler Einsatz

Curetas, subgingival

Curettes, subgingival

Curette, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia



Langer	Langer	Langer
SP29-142-02	SP29-142-04	SP29-142-06
SP29-143-02	SP29-143-04	SP29-143-06

Fig. 1/2



Mandible, lateral teeth,
mesial / distal




Fig. 3/4



Maxilla, lateral teeth,
mesial / distal



Fig. 5/6



Maxilla and mandible,
incisors, universal

Langer

SP29-142-02 -

SP29-143-06

17,5 cm / 6 7/8"

Universal curette

hollow handle, 8 mm

GR = grey

Y = yellow

R = red

G = green

V = violet

W = white

B = blue

N = black




Roncati Scaler
SP29-144-00
 17,5 cm / 6 7/8"
 Scaler
 hollow handle, 10 mm



Roncati
SP29-144-01
 17,5 cm / 6 7/8"
 Curette, micro, straight
 hollow handle, 10 mm



Roncati
SP29-144-02
 17,5 cm / 6 7/8"
 Curette, micro, curved
 hollow handle, 10 mm



Roncati
SP29-145-00
 17,5 cm / 6 7/8"
 Curette, micro, curved
 hollow handle, 10 mm

Ti

Scaling- Instruments

Scaling Instruments
Instruments de Détartrage
Scaling-Instrumente

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia

Why scaling instruments should always be sharp:


1. Sharp instruments allow precision work with little application of force.
2. Blunt instruments polish the dental stone, bruise the soft tissue and require more application of force.

What you need for instrument sharpening:

1. Sharpening stones (rough, fine, and round for burr removal)
2. Sharpening oil (prevents excessive heat build-up; preserves the sharpening capacity of the sharpening stone)
3. Test stick (made of acrylic, for sharpness test)
4. Good light (sharp edges do not reflect light)
5. Clamping device (e.g. vise)

What to check after sharpening:

1. Roundness of back surfaces
2. Continuous sharpness of cutting edges
3. Perfect shape of working end (e.g. check against catalog illustration)
4. After several sharpening cycles, verify that the instrument is still strong enough for application. Replace instrument if necessary.



Color coding of scalers and Gracey curettes

1. All scalers are color-coded green – 3 points distributed across the diameter. At the same time, the color code marks the finer working end.
2. The Gracey Standard, Gracey Rigid and Gracey Mini curettes also have the color code on the side of the finer working end.
Example: Gracey 1/2 = color code applied on side of working end no. 1.
3. For even better differentiation between Gracey Standard and Gracey Mini on the one hand and Gracey Rigid on the other hand, Gracey Rigid curettes are coded with six color points.
4. Color-coded Gracey Standard curettes are also available with a 10 mm handle.

left: color
(lower figure,
e.g. insert 3)



right: no color
(higher figure,
e.g. insert 4)



	Ø		
Gracey Mini	8 mm	x	
Gracey Standard	8 mm	x	
Gracey Big	10 mm	x	
Gracey Rigid	8 mm		x



Fields of application for Gracey curettes

-  1-2,3-4, 5-6 Incisors
-  7-8, 9-10 Premolars and molars (buccal/lingual)
-  11-12, 15-16 Premolars and molars, mesial Furcation: mesial portion of the distal root
-  13-14, 17-18 Premolars and molars, distal Furcation: distal portion of the mesial root

Gracey mini curettes, subgingival

Gracey Mini-Küretten, subgingivaler Einsatz

Gracey mini curetas, subgingival

Gracey mini curettes, subgingival

Gracey mini curette, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-146-02

Fig. 1/2



Incisors, universal,
first shaft longer, working
end shorter




SP29-146-04

Fig. 3/4



Incisors, premolars,
universal,
first shaft longer,
working end shorter




SP29-146-06

Fig. 5/6




Incisors, premolars,
universal,
first shaft longer,
working end shorter



SP29-146-08

Fig. 7/8



Premolars, molars,
lingual/buccal,
first shaft longer, working
end shorter

Gracey Mini

SP29-146-02 -

SP29-146-18

17,5 cm / 6 7/8"

hollow handle, 8 mm

GR = grey

Y = yellow

R = red


G = green

V = violet

W = white


B = blue

N = black




SP29-146-10

Fig. 9/10




Premolars, molars,
lingual/buccal,
first shaft longer,
working end shorter




SP29-146-12

Fig. 11/12



Premolars, molars,
mesial, mesiolingual,
mesiobuccal,
first shaft longer,
working end shorter




SP29-146-14

Fig. 13/14




Premolars, molars, distal,
first shaft longer,
working end shorter




SP29-146-16

Fig. 15/16



Premolars, molars, mesial,
first shaft longer,
working end shorter



SP29-146-18

Fig. 17/18



Premolars, molars, distal,
first shaft longer,
working end shorter

Gracey standard curettes, subgingival

Gracey Standard-Küretten, subgingivaler Einsatz

Gracey curetas estándar, subgingival

Gracey curettes standard, subgingival

Gracey curette standard, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología


Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-147-02

Fig. 1/2




Incisors, universal




SP29-147-04

Fig. 3/4




Incisors, premolars,
universal




SP29-147-06

Fig. 5/6



Incisors, premolars,
universal



SP29-147-08

Fig. 7/8



Premolars, molars,
lingual / buccal

Gracey Standard


SP29-147-02 -

SP29-147-18


17,5 cm / 6 7/8"

hollow handle, 8 mm






SP29-147-10
Fig. 9/10




Premolars, molars,
lingual / buccal



SP29-147-12
Fig. 11/12




Premolars, molars, mesial,
mesiobuccal



SP29-147-14
Fig. 13/14




Premolars, molars, distal




SP29-147-16
Fig. 15/16



Premolars, molars, mesial



SP29-147-18
Fig. 17/18



Premolars, molars, distal

Gracey Rigid curettes, subgingival

Gracey Rigid-Küretten, subgingivaler Einsatz

Gracey Rigid curetas, subgingival

Gracey Rigid curettes, subgingival

Gracey Rigid curette, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-148-02

Fig. 1/2



Incisors, universal,
extra rigid, for stubborn
dental plaque




SP29-148-04

Fig. 3/4



Incisors, premolars,
universal,
extra rigid, for stubborn
dental plaque




SP29-148-06

Fig. 5/6



Incisors, premolars,
universal,
extra rigid, for stubborn
dental plaque



SP29-148-08

Fig. 7/8




Premolars, molars, lingual,
buccal,
extra rigid, for stubborn
dental plaque

Gracey Rigid SP29-148-02 - SP29-148-18

17,5 cm / 6 7/8"


hollow handle, 8 mm






SP29-148-10

Fig. 9/10




SP29-148-12

Fig. 11/12




SP29-148-14

Fig. 13/14




SP29-148-16

Fig. 15/16



SP29-148-18

Fig. 17/18



Premolars, molars, lingual,
 buccal,
 extra rigid, for stubborn
 dental plaque




Premolars, molars, mesial,
 mesiolingual, mesiobuccal,
 extra rigid, for stubborn
 dental plaque



Premolars, molars,
 distal
 extra rigid, for stubborn
 dental plaque



Premolars, molars, mesial,
 extra rigid, for stubborn
 dental plaque




Premolars, molars,
 distal
 extra rigid, for stubborn
 dental plaque

Gracey curettes with 10 mm handle

Gracey Küretten mit 10 mm-Griff
 Gracey curetas para mangos 10 mm
 Gracey curettes avec manches à 10 mm
 Gracey curette con manici 10 mm

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
 Limpiadores de sarro, parodontología
 Instruments à détartrer, parodontologie
 Ablatori di tartaro, parodontologia




SP29-149-02

Fig. 1/2



Incisors, universal




SP29-149-04

Fig. 3/4



Incisors, premolars, universal




SP29-149-06

Fig. 5/6



Incisors, premolars, universal



SP29-149-08

Fig. 7/8



Premolars, molars, lingual, buccal


Big Gracey

**SP29-149-02 -
 SP29-149-18**

17,5 cm / 6 7/8"


hollow handle, 10 mm






SP29-149-10

Fig. 9/10




SP29-149-12

Fig. 11/12




SP29-149-14

Fig. 13/14




SP29-149-16

Fig. 15/16



SP29-149-18


Fig. 17/18



Premolars, molars, lingual,
buccal




Premolars, molars,
mesial, mesiolingual,
mesiobuccal,



Premolars, molars,
distal



Premolars, molars, mesial



Premolars, molars,
distal

Curettes, subgingival

Küretten, subgingivaler Einsatz
Curetas, subgingival
Curettes, subgingival
Curette, subgingivale

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia



SP29-150-01

Fig. 1



Furcation curette
0,9 mm
buccal / lingual



SP29-150-02

Fig. 2



Furcation curette
3 mm
buccal / lingual



SP29-150-03

Fig. 3



Furcation curette
0,9 mm
mesial / distal



SP29-150-04


Fig. 4



Furcation curette
1,3 mm
mesial / distal

**SP29-150-01 -
SP29-150-04**

17,5 cm / 6 7/8"
hollow handle, 8 mm




SP29-150-05

Fig. 5



Furcation curette
 0,9 mm
 labial / lingual



Prichard
SP29-151-01

Fig. 1/2



Bone curette,
 extra long, robust
 construction, double
 cutting

Gracey mini curettes, subgingival

Gracey Mini-Küretten, subgingivaler Einsatz

Gracey mini curetas, subgingival

Gracey mini curettes, subgingival

Gracey mini curette, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-152-06

Ti

Fig. 5/6

17,5 cm / 6 $\frac{7}{8}$ "

hollow handle, 10 mm




SP29-152-08

Ti

Fig. 7/8


17,5 cm / 6 $\frac{7}{8}$ "

hollow handle, 10 mm



Premolars, molars,
mesial, mesiolingual,
mesiobuccal,
first shaft longer, working
end shorter

Premolari, molari,
mesiali, mesiolinguali,
mesiobuccali,
1° stelo più lungo, estremità
di lavoro più breve




SP29-152-12

Ti


Fig. 11/12

17,5 cm / 6 7/8"

hollow handle, 10 mm



Premolars, molars, distal,
first shaft longer, working
end shorter



SP29-152-14

Ti

Fig. 13/14

17,5 cm / 6 7/8"

hollow handle, 10 mm

Gracey standard curettes, subgingival

Gracey Standard-Küretten, subgingivaler Einsatz

Gracey curetas estándar, subgingival

Gracey curettes standard, subgingival

Gracey curette standard, subgingivale


Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie

Ablatori di tartaro, parodontologia




SP29-153-06

Ti

Fig. 5/6

17,5 cm / 6 7/8"

hollow handle, 10 mm




SP29-153-08

Ti

Fig. 7/8

17,5 cm / 6 7/8"

hollow handle, 10 mm




SP29-153-12

Ti

Fig. 11/12

17,5 cm / 6 7/8"

hollow handle, 10 mm



SP29-153-14

Ti

Fig. 13/14

17,5 cm / 6 7/8"

hollow handle, 10 mm

Sharpening stones

Schleifsteine
Piedras de afilar
Meules
Pietre per affilare

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia



Record Arkansas SP29-154-00

115 x 40 x 10/2 mm
multi shape, white

For pregrinding of blunt instruments. The rounded edges and bevelled surfaces enable the sharpening of radius and surfaces



Original Arkansas SP29-154-01

100 x 25 x 7 mm
flat, grey-white

For repeated grinding and finishing of scalers and curettes



Original Arkansas SP29-154-02

100 mm lang, Ø 8 mm
round, grey-white

For burr removing from cutting edges



India-Korund SP29-154-05

115 x 45 x 2 - 10 mm

delicate, brown-beige
Para preafilar y volverles a dar forma a instrumentos



India-Korund SP29-154-06

100 mm lang, Ø 5 - 10 mm

round, delicate, brown-beige
For burr removing from cutting edges



SP29-154-08

75 mm lang, Ø 6,0 mm

Test stick plastic
package of 6 pieces
For testing the sharpness of scalers and curettes



Record Arkansas

SP29-155-00

150 x 50 mm
in wooden case



SP29-154-09

Sharpening stone oil

Periodontal probes

Parodontometer
Sondas parodontales
Sondes parodontales
Sonde parodontali

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia




SP29-156-01

Iglhaut
SP29-156-01
17,5 cm / 6 7/8"
solid handle, 6 mm



SP29-157-40

Nabers
SP29-157-40
17,5 cm / 6 7/8"
hollow handle, 8 mm
Furcation probe, graduated working end



SP29-158-08

Fig. 8
Grad. 3-6-8-11

SP29-158-10

Fig. 10
Grad. 1-2-3-5-7-8-9-10

SP29-158-11


Fig. 11
Grad. 3-6-8-11

SP29-159-00

Fig. W
Grad. 1-2-3-5-7-8-9-10

SP29-158-08 - SP29-159-00

16 cm / 6 1/8"
solid handle, 6 mm



SP29-160-11	SP29-160-13	SP29-160-21	SP29-160-23
SP29-161-11	SP29-161-13	SP29-161-21	SP29-161-23
Fig. CPG 11,5 WHO Grad. 3,5-5,5-8,5-11,5 with ball	Fig. CPG 12 Grad. 3-6-9-12	Fig. CP 15 UNC Grad. 1-2-3-4-...-14-15	Fig. CPNG 22 Grad. 2-4-6-8-10-12




Probe handles see page 244

Periodontal probes

Parodontometer
Sondas parodontales
Sondes parodontales
Sonde parodontali

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia





0,8 mm

1/1

Williams-Goldman-Fox

SP29-163-02

Fig. 1-2-3-5-7-8-9-10

0,6 mm

SP29-163-04

Fig. CP 11/23

Grad. 3-6-8-11

0,6 mm

SP29-163-06

Fig. CP 12/23

Grad. 3-6-9-12

0,5 mm

0,8 mm

SP29-163-08

Fig. GF/3

Grad. 1-2-3-5-7-8-9-10





0,5 mm

0,8 mm

SP29-163-10

Fig. UM/3

Grad. 3-6-8



0,5 mm

0,8 mm

1/1

SP29-163-12

Fig. W/3

Grad. 1-2-3-5-7-8-9-10

0,6 mm

SP29-163-14

Fig. GF/6

Grad. 1-2-3-5-7-8-9-10

0,6 mm

SP29-163-16

Fig. UM/6

Grad. 3-6-8

0,5 mm

SP29-163-18

Fig. W/6






Grad. 1-2-3-5-7-8-9-10

0,6 mm

SP29-163-20

Fig. GF/17

Grad. 1-2-3-5-7-8-9-10



SP29-163-02 -

SP29-163-20

17,5 cm / 6 7/8"

hollow handle, 8 mm

Periodontal probes

Parodontometer
Sondas parodontales
Sondes parodontales
Sonde parodontali

Scalers, periodontology


Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia



SP29-163-22

Fig. UM/17


Grad. 3-6-8



SP29-163-24

Fig. W/17


Grad. 1-2-3-5-7-8-9-10



SP29-163-26

Fig. GF/23



Grad. 1-2-3-5-7-8-9-10



SP29-163-28

Fig. UM/23


Grad. 3-6-8



SP29-163-30

Fig. W/23

Grad. 1-2-3-5-7-8-9-10




SP29-163-22 - SP29-163-30

17,5 cm / 6 7/8"

hollow handle, 8 mm




Cross
SP29-164-10
17,5 cm / 6 7/8"
hollow handle, 8 mm



Zürich
SP29-165-02
17,5 cm / 6 7/8"
hollow handle, 8 mm




FinoLine



SP29-165-03
Fig. 3




SP29-165-06
Fig. 6



SP29-165-08
Fig. 8 A




SP29-165-09
Fig. 9




SP29-165-12
Fig. 12




SP29-165-17
Fig. 17




SP29-165-23
Fig. 23



SP29-165-30
Fig. 3A




SP29-165-60
Fig. 6XL



White
SP29-165-12w
Fig. 12

SP29-165-03 -
SP29-165-12w
17,5 cm / 6 7/8"
extra delicate tooth probe
solid handle, 6 mm



SP29-166-03

Fig. 3


Detection of furcations and for easy entry into interproximal spaces



SP29-166-06

Fig. 6

Detection of plaque, caries, root canals and furcations



SP29-166-08

Fig. 8A




SP29-166-09

Fig. 9



Briault
SP29-166-11

Fig. EXD 11



SP29-166-12


Fig. EXD 11-12

SP29-166-03 -
SP29-166-12
17,5 cm / 6 7/8"
solid handle, 6 mm




White
SP29-166-17
Fig. 17

Mainly for sub-gingival examination and approx. caries




White
SP29-166-23
Fig. 23


For general examination, especially for incisors



SP29-166-33
Fig. 3A



SP29-166-54
Fig. 54




SP29-166-66
Fig. 6XL

Detection of plaque, caries, root canals and furcations

SP29-166-17 -
SP29-166-66
17,5 cm / 6 7/8"
solid handle, 6 mm




**SP29-167-01 -
SP29-167-23**
17,5 cm / 6 7/8"
hollow handle, 8 mm



SP29-167-01
Fig. UT 1/2




SP29-167-02
Fig. UT 2/2



SP29-167-03
Fig. 3CH/3CH



SP29-167-04
Fig. 6/3




SP29-167-06
Fig. 23/6







SP29-167-16
Fig. 16/17




SP29-167-17
Fig. 17/23




SP29-167-23
Fig. 23/17A






Nabers
SP29-168-01
Fig. 1






Nabers
SP29-168-02
Fig. 2




SP29-168-05
Fig. 5/6



SP29-168-23
Fig. 16/23



SP29-168-01 -
SP29-168-23
17,5 cm / 6 7/8"
hollow handle, 8 mm



SP29-169-08

Fig. EXG 8

SP29-169-09

Fig. EXG 9

SP29-169-10

Fig. EXG 10

SP29-169-23


Fig. EXG 23




Probe handles



**SP29-170-01 -
 SP29-171-34**
 17,5 cm / 6 7/8"
 hollow handle, 8 mm
 double ended




SP29-170-01
 Fig. GM 1A




SP29-170-02
 Fig. GM 2A



SP29-171-00
 Fig. GM 0



SP29-171-01
 Fig. GM 1



grey



yellow



red



green



violet




white



blue




light-blue



SP29-171-02

Fig. GM 2




SP29-171-03

Fig. GM 3



SP29-171-04

Fig. GM 4



Zürich
SP29-171-34

Ø 1,9 mm





**SP29-172-00 -
 SP29-172-44**
 17,5 cm / 6 7/8"
 hollow handle, 8 mm
 double ended



SP29-172-00
 Fig. GM 125/126



SP29-172-02
 Fig. GM 127/128



SP29-172-04
 Fig. GM 129/130



SP29-172-06
 Fig. GM 133/134



SP29-172-08
 Fig. GM 153/154



SP29-172-10
 Fig. GM 14



SP29-172-12
 Fig. GM 38/39



SP29-172-14
 Fig. GM E1



SP29-172-16
 Fig. GM E2





SP29-172-18
 Fig. GM E3



SP29-172-20
 Fig. GM 6



SP29-172-22
 Fig. GM 17S






SP29-172-24
Fig. GM 17



SP29-172-26
Fig. GM 18



SP29-172-28
Fig. GM 19



SP29-172-30
Fig. GM 17W



SP29-172-32
Fig. GM 18W



SP29-172-34
Fig. GM 19W




SP29-172-35
Fig. GM 61/62




SP29-172-36
Fig. GM 63/64




SP29-172-38
Fig. GM 65/66




SP29-172-40
Fig. GM 31L



SP29-172-42
Fig. GM 32L



SP29-172-44
Fig. GM 33L



Enamel chisels

Schmelzmeißel
Cinceles para esmalte
Ciseaux à émail
Tagliasmalto

Scalers, periodontology

Reinigungsinstrumente, Parodontologie
Limpiadores de sarro, parodontología
Instruments à détartrer, parodontologie
Ablatori di tartaro, parodontologia




1/1

Black


**SP29-173-01 -
SP29-173-80**

17,5 cm / 6 7/8"
solid handle, 6 mm
double ended




SP29-173-01

Fig. 1




SP29-173-02

Fig. 2



SP29-173-26

Fig. 26



SP29-173-27

Fig. 27



SP29-173-28

Fig. 28



SP29-173-29

Fig. 29



SP29-173-52

Fig. 51/52



SP29-173-54

Fig. 53/54



SP29-173-64

Fig. 63/64



SP29-173-66

Fig. 65/66



SP29-173-78

Fig. 77/78




SP29-173-80

Fig. 79/80




1/1








Wedelstaedt
SP29-174-02 -
SP29-174-06
17,5 cm / 6 7/8"




SP29-174-02
Fig. 1/2



SP29-174-04
Fig. 3/4



SP29-174-06
Fig. 5/6



solid handle, 6 mm
double ended

Root canal instruments

Nervkanal-Instrumente
Instrumentos para canales

Instrumenti a canaux

Strumenti per trattamenti canalari

Scalers, periodontology

Reinigungsinstrumente, Parodontologie

Limpiadores de sarro, parodontología

Instrumenti a détartrer, parodontologie

Ablatori di tartaro, parodontologia



1/1


SP29-175-31 -

SP29-176-01

17,5 cm / 6 7/8"

hollow handle, 8 mm

double ended




Zürich

SP29-175-31

Fig. 30/40



Pluggen




Zürich

SP29-175-51

Fig. 50/90



Pluggen



Stewart

SP29-176-01

Endo probe



grey



yellow



red



green



violet




white



blue




black




SP29-177-03

Fig. 30




SP29-177-04

Fig. 40




SP29-177-05

Fig. 50




SP29-177-06

Fig. 60




SP29-177-09

Fig. 90



SP29-177-10

Fig. 100



SP29-177-12

Fig. 120

Plugger



1/4

SP29-177-03 -

SP29-177-12

17,5 cm / 6 7/8"
solid handle, 6 mm
single ended

Root canal instruments

Nervkanal-Instrumente

Instrumentos para canales

Instrumenti à canaux

Strumenti per trattamenti canalari


Scalers, periodontology

Reinigungsinstrumente, Parodontologie


Limpiadores de sarro, parodontología

Instrumenti à détartrer, parodontologie


Ablatori di tartaro, parodontologia



GR grey Y yellow R red G green V violet W white B blue N black



SP29-180-37
Fig. 7



SP29-180-38
Fig. 8




SP29-180-39
Fig. 9



SP29-180-40
Fig. 10



SP29-180-41
Fig. 11



**SP29-180-37 -
 SP29-180-41**
 17,5 cm / 6 7/8"
 solid handle, 6 mm
 double ended

Gingivectomy knives

Gingivektomiemesser

Bisturís para la gingivectomía

Bistouris pour la gingivectomie

Bisturie per gingivectomia


Scalers, periodontology

Reinigungsinstrumente, Parodontologie



Limpiadores de sarro, parodontología

Instruments à détartrer, parodontologie


Ablatori di tartaro, parodontologia



SP29-181-03



SP29-181-04



Allen-Iglhaut


SP29-181-03 -

SP29-181-04

16 cm / 6 4/8"


hollow handle, 10 mm

double ended




Kirkland
SP29-182-51


For gingivectomy or gingivoplastic surge
Very effective in retro molar region



Kirkland
SP29-182-61




Kirkland
SP29-183-61




Orban
SP29-184-01

Fig. 1



Orban
SP29-184-02

Fig. 2



Orban
SP29-184-12

Fig. 1/2



SP29-182-51 -
SP29-184-12
16 cm / 6 1/8"
solid handle, 6 mm
single ended

Counter angled shank, especially suitable for excision of
interproximal tissue in gingivectomy

Integrity in Every Incision.

Thank you for taking the time to learn about Surgiplex SRL – your trusted partner in surgical innovation and healthcare supply.

Head Office

Strada Vanatori 3, Pravalia 6,
Sector 5, Bucharest, Romania

Email: surgiplexsrl@gmail.com

WhatsApp: +92 333 8635932

EU Contact: +40 765 651 717

GCC Contact: +974 3397 9907

Website: www.surgiplexsrl.com

SURGIPILEX